

LINEAR SHAKER

Model: 264



FEATURES

- The Scientific linear shakers create a shaking motion ideal for mixing liquids in flasks and conducive to culturing cells and for testing cans/bottles in the beverage industry
- The robust, functional design ensures reliable trouble free operation
- Microprocessor based speed control with easy to use RPM setting and digital indication of speed
- RPM: 50-400
- Digital timer setting with the option of run-down or continuous operation
- Loose rubber mat fitted to the platform
- 20kg load capacity

NOTE

- A full range of accessories are available for all clamping requirements
- Unit supplied complete with platform and 4 clamping bars
- Optional beverage platform available



ALL SCIENTIFIC EQUIPMENT COMES WITH A TWO
YEAR WARRANTY ON COMPONENTS AND
DEFECTIVE WORKMANSHIP - Ts & Cs APPLY.

www.scientific.co.za

LINEAR SHAKER


Model: 264

SPECIFICATIONS

POWER	230W
VOLTAGE	230V / 50Hz
PACKAGING (W x D x H)	675 x 600 x 265 (mm)
UNIT (W x D x H)	650 x 570 x 250 (mm)
NET WEIGHT	20kg
GROSS WEIGHT	23kg
SHIPPING VOLUME	0.11m ³
PLATFORM DIMENSIONS	570 x 460 (mm)
STROKE	40mm
RPM	50 - 400
LOAD CAPACITY	20kg

OPTIONAL EXTRAS

CODE	DESCRIPTION
269A	50ml Flask Clamp (Capacity 30)
269B	125ml Flask Clamp (Capacity 30)
269C	250ml Flask Clamp (Capacity 20)
269D	500ml Flask Clamp (Capacity 12)
269E	1000ml Flask Clamp (Capacity 4)
269F	2000ml Flask Clamp (Capacity 2)



NOTE

- When viewing the appliance from the front in its usual operating position, the width (W) of the product is the total distance from left to right; the depth (D) of the product is the total distance from the front to the back; the height (H) is the total distance from the bottom of the product to the top. Pictures may vary from actual product.

DISCLAIMER

- Output performance figures quoted are dependant on various factors. Scientific reserves the right, without notice, to make changes and revisions to product specifications, materials and design, when we believe it will provide better performance, durability, and efficiency. Unit measurements may vary depending on point of reference.

WARRANTY

- All Scientific equipment comes with a standard ONE year warranty, The 1 Year Extended Warranty is only valid for customers who have registered their warranty online at www.scientific.co.za within the first year after making their purchase.