



## HOTPLATES

Model: 501



## FEATURES

- Laboratory hotplates are ideally suited for applications where high temperature stability is not critical
- Heavy duty mild steel heating top with special Teflon finish
- Mild steel body treated for long term resistance to corrosion with an epoxy powder coated finish
- Unique body design to ensure maximum cooling airflow to the controls and wiring
- Suitable for continuous operation
- Variable energy regulator for optimum heat control
- Pilot light for power indication
- Maximum operating temperature: 250°C

### NOTE

- Sandbath attachments available
- Optional remote control available to isolate control from the unit to the outside of a fumehood



ALL SCIENTIFIC EQUIPMENT COMES WITH A TWO  
YEAR WARRANTY ON COMPONENTS AND  
DEFECTIVE WORKMANSHIP - Ts & Cs APPLY.

[www.scientific.co.za](http://www.scientific.co.za)



## HOTPLATES



Model: 501

## SPECIFICATIONS

POWER	2000W
VOLTAGE	230V / 50Hz
PACKAGING (W x D x H)	550 x 400 x 255 (mm)
UNIT (W x D x H)	455 x 340 x 165 (mm)
NET WEIGHT	22kg
GROSS WEIGHT	23kg
SHIPPING VOLUME	0.06m <sup>3</sup>
SURFACE DIMENSIONS	450 x 300 (mm)
MAX TEMPERATURE	250°C

## OPTIONAL EXTRAS

CODE	DESCRIPTION
510	Sandbath Attachment
501E	Remote control



### NOTE

- When viewing the appliance from the front in its usual operating position, the width (W) of the product is the total distance from left to right; the depth (D) of the product is the total distance from the front to the back; the height (H) is the total distance from the bottom of the product to the top. Pictures may vary from actual product.

### DISCLAIMER

- Output performance figures quoted are dependent on various factors. Scientific reserves the right, without notice, to make changes and revisions to product specifications, materials and design, when we believe it will provide better performance, durability, and efficiency. Unit measurements may vary depending on point of reference.

### WARRANTY

- All Scientific equipment comes with a standard ONE year warranty. The 1 Year Extended Warranty is only valid for customers who have registered their warranty online at [www.scientific.co.za](http://www.scientific.co.za) within the first year after making their purchase.