



LABORATORY TROLLEY

Model: 552



FEATURES

- Laboratory and general purpose trolleys are characterized by their robust stainless steel construction and hygienic, easy to clean design
- Wheels, with vibration absorbing tyres, ensure ease of manoeuvrability and no snagging on floor joints or lift entrances
- The top tray is at a convenient working height to allow for easy transfer of items to and from any standard bench
- Units are easy to assemble, instructions provided
- 100mm easy glide, silent running, heavy duty castors
- Carrying capacity 50 kg per shelf
- Sturdy stainless steel framework

NOTE

- Supplied in knock down form
- Supplied with 2 stainless steel trays
- Non marking tyres to ensure floors remain clean



ALL SCIENTIFIC EQUIPMENT COMES WITH A TWO
YEAR WARRANTY ON COMPONENTS AND
DEFECTIVE WORKMANSHIP - Ts & Cs APPLY.

www.scientific.co.za



LABORATORY TROLLEY



Model: 552

SPECIFICATIONS

PACKAGING (W x D x H)	1000 x 205 x 600 (mm)
UNIT (W x D x H)	930 x 550 x 930 (mm)
NET WEIGHT	20kg
GROSS WEIGHT	21kg
SHIPPING VOLUME	0.12m ³
NUMBER OF TRAYS	2
DISTANCE BETWEEN TRAYS	450mm

OPTIONAL EXTRAS

CODE	DESCRIPTION
561	Extra stainless steel tray



NOTE

- When viewing the appliance from the front in its usual operating position, the width (W) of the product is the total distance from left to right; the depth (D) of the product is the total distance from the front to the back; the height (H) is the total distance from the bottom of the product to the top. Pictures may vary from actual product.

DISCLAIMER

- Output performance figures quoted are dependent on various factors. Scientific reserves the right, without notice, to make changes and revisions to product specifications, materials and design, when we believe it will provide better performance, durability, and efficiency. Unit measurements may vary depending on point of reference.

WARRANTY

- All Scientific equipment comes with a standard ONE year warranty. The 1 Year Extended Warranty is only valid for customers who have registered their warranty online at www.scientific.co.za within the first year after making their purchase.



ALL SCIENTIFIC EQUIPMENT COMES WITH A TWO
YEAR WARRANTY ON COMPONENTS AND
DEFECTIVE WORKMANSHIP - Ts & Cs APPLY.

www.scientific.co.za

Published Date: 2020.09.25