



## 1.5L BLENDER

Model: 722



## FEATURES

- Scientific blenders are durable and robust for use in agricultural, microbiological and cement laboratories
- Multi-function
- Two-in-one stainless steel cutting blade
- Industrial high speed motor
- Premium quality 1.5L polycarbonate jug with mixing rod
- Precision technology



ALL SCIENTIFIC EQUIPMENT COMES WITH A TWO YEAR WARRANTY ON COMPONENTS AND DEFECTIVE WORKMANSHIP - Ts & Cs APPLY.

[www.scientific.co.za](http://www.scientific.co.za)






## 1.5L BLENDER

Model: 722

### SPECIFICATIONS

POWER	1200W
VOLTAGE	230V / 50Hz
PACKAGING (W x D x H)	410 x 215 x 240 (mm)
UNIT (W x D x H)	200 x 225 x 420 (mm)
NET WEIGHT	4.5kg
GROSS WEIGHT	5.5kg
SHIPPING VOLUME	0.04m <sup>3</sup>
JUG CAPACITY	1.5ℓ



#### NOTE

- When viewing the appliance from the front in its usual operating position, the width (W) of the product is the total distance from left to right; the depth (D) of the product is the total distance from the front to the back; the height (H) is the total distance from the bottom of the product to the top. Pictures may vary from actual product.

#### DISCLAIMER

- Output performance figures quoted are dependent on various factors. Scientific reserves the right, without notice, to make changes and revisions to product specifications, materials and design, when we believe it will provide better performance, durability, and efficiency. Unit measurements may vary depending on point of reference.

#### WARRANTY

- All Scientific equipment comes with a standard ONE year warranty, The 1 Year Extended Warranty is only valid for customers who have registered their warranty online at [www.scientific.co.za](http://www.scientific.co.za) within the first year after making their purchase.