

600mm LABORATORY HOTPLATES

Model: 732



FEATURES

- ❑ Ideally suited for applications where even temperature distribution is required
- ❑ Heavy duty aluminium top
- ❑ Excellent heat distribution
- ❑ Digital heat control for accurate temperature setting
- ❑ Maximum temperature setting of 350°C




ALL SCIENTECH EQUIPMENT COMES WITH A ONE
YEAR WARRANTY ON COMPONENTS AND
DEFECTIVE WORKMANSHIP

600mm LABORATORY HOTPLATES

Model: 732

SPECIFICATIONS

POWER	3000W
VOLTAGE	230V / 50Hz
PACKAGING (W x D x H)	675 x 620 x 245 (mm)
UNIT (W x D x H)	600 x 550 x 150 (mm)
NET WEIGHT	31kg
GROSS WEIGHT	34kg
SHIPPING VOLUME	0.1m ³
MAX TEMPERATURE	350°C
PLATE SIZE	600 x 450 (mm)



TERMS & CONDITIONS

NOTE: Pictures may vary from actual product. | **DISCLAIMER:** Output performance figures quoted are dependent on various factors. Scientific reserves the right, without notice, to make changes and revisions to product specifications, materials and design, when we believe it will provide better performance, durability, and efficiency. Unit measurements may vary depending on point of reference. | **WARRANTY:** All Scientech equipment comes with a standard ONE year warranty.

REMEMBER TO COMPLETE THE ONLINE WARRANTY REGISTRATION



ALL SCIENTECH EQUIPMENT COMES WITH A ONE
YEAR WARRANTY ON COMPONENTS AND
DEFECTIVE WORKMANSHIP

Published Date: 2020.09.25