



GLASSWARE DRIER

Model: 310



FEATURES

- This simple, robust unit is ideally suited for the rapid and efficient drying of laboratory glassware
- Grade 304 stainless steel body construction
- Reliable incoloy heating elements
- Toughened glass sliding doors on both sides for ease of access
- Variable heat control via energy regulator

NOTE

- Unit provided with 2 chrome plated wire shelves



ALL SCIENTIFIC EQUIPMENT COMES WITH A TWO
YEAR WARRANTY ON COMPONENTS AND
DEFECTIVE WORKMANSHIP - Ts & Cs APPLY.

www.scientific.co.za




GLASSWARE DRIER

Model: 310

SPECIFICATIONS

POWER	1000W
VOLTAGE	230V / 50Hz
PACKAGING (W x D x H)	1120 x 620 x 585 (mm)
UNIT (W x D x H)	1040 x 555 x 530 (mm)
NET WEIGHT	27kg
GROSS WEIGHT	30kg
SHIPPING VOLUME	0.41m ³
CHAMBER DIMENSIONS	965 x 495 x 415 (mm)
CHAMBER VOLUME	180ℓ



NOTE

- When viewing the appliance from the front in its usual operating position, the width (W) of the product is the total distance from left to right; the depth (D) of the product is the total distance from the front to the back; the height (H) is the total distance from the bottom of the product to the top. Pictures may vary from actual product.

DISCLAIMER

- Output performance figures quoted are dependent on various factors. Scientific reserves the right, without notice, to make changes and revisions to product specifications, materials and design, when we believe it will provide better performance, durability, and efficiency. Unit measurements may vary depending on point of reference.

WARRANTY

- All Scientific equipment comes with a standard ONE year warranty. The 1 Year Extended Warranty is only valid for customers who have registered their warranty online at www.scientific.co.za within the first year after making their purchase.



ALL SCIENTIFIC EQUIPMENT COMES WITH A TWO
YEAR WARRANTY ON COMPONENTS AND
DEFECTIVE WORKMANSHIP - Ts & Cs APPLY.

www.scientific.co.za

Published Date: 2026.09.16