


40L DESICCATOR CABINET

Model: 710



FEATURES

- Scientific laboratory desiccators are sealable enclosures containing desiccants used for preserving moisture-sensitive items
- For indoor use only
- Stainless steel 304 interior
- Epoxy powder coated exterior
- Magnetic door closure with toughened glass window

NOTE

- Supplied with two perforated trays and one tray for desiccant
- Desiccants are not included



ALL SCIENTIFIC EQUIPMENT COMES WITH A TWO
YEAR WARRANTY ON COMPONENTS AND
DEFECTIVE WORKMANSHIP - Ts & Cs APPLY.

www.scientific.co.za




40L DESICCATOR CABINET

Model: 710

SPECIFICATIONS

PACKAGING (W x D x H)	590 x 575 x 550 (mm)
UNIT (W x D x H)	480 x 480 x 465 (mm)
NET WEIGHT	8kg
GROSS WEIGHT	9kg
SHIPPING VOLUME	0.19m ³
CHAMBER DIMENSIONS	350 x 330 x 350 (mm)
CHAMBER VOLUME	40ℓ



NOTE

- When viewing the appliance from the front in its usual operating position, the width (W) of the product is the total distance from left to right; the depth (D) of the product is the total distance from the front to the back; the height (H) is the total distance from the bottom of the product to the top. Pictures may vary from actual product.

DISCLAIMER

- Output performance figures quoted are dependent on various factors. Scientific reserves the right, without notice, to make changes and revisions to product specifications, materials and design, when we believe it will provide better performance, durability, and efficiency. Unit measurements may vary depending on point of reference.

WARRANTY

- All Scientific equipment comes with a standard ONE year warranty. The 1 Year Extended Warranty is only valid for customers who have registered their warranty online at www.scientific.co.za within the first year after making their purchase.



**ALL SCIENTIFIC EQUIPMENT COMES WITH A TWO
YEAR WARRANTY ON COMPONENTS AND
DEFECTIVE WORKMANSHIP - Ts & Cs APPLY.**

www.scientific.co.za

Published Date: 2026.09.17