
Scientific | Laboratory Equipment


PROUDLY MANUFACTURED IN SOUTH AFRICA

MAGNETIC STIRRER
Model: 105









FEATURES

- Magnetic Stirrers are suitable for most types of laboratory applications where controlled stirring and heating is required
- Aluminum top for excellent heat distribution
- Reliable heating elements
- Temperature range 50°C - 350°C, controlled by analogue potentiometer
- Reliable solid state electronic speed control
- Stirring bar 40 x 8mm
- RPM: 100 - 1800

NOTE

- Supplied with probe stand, clamp and ceramic stirrer bar

QUALITY GUARANTEED

ALL SCIENTIFIC EQUIPMENT COMES WITH A TWO YEAR WARRANTY ON COMPONENTS AND DEFECTIVE WORKMANSHIP - Ts & Cs APPLY.


Scientific | Laboratory Equipment

PROUDLY MANUFACTURED IN SOUTH AFRICA

MAGNETIC STIRRER
Model: 105

SPECIFICATIONS

POWER	500W
VOLTAGE	230V / 50Hz
PACKAGING (W x D x H)	230 x 310 x 215 (mm)
UNIT (W x D x H)	180 x 270 x 330 (mm)
NET WEIGHT	4.6kg
GROSS WEIGHT	5kg
SHIPPING VOLUME	0.02m^3
LOAD CAPACITY	5ℓ
RPM	100 - 1800
PLATE SIZE	170 x 170 (mm)



NOTE
- When viewing the appliance from the front in it's usual operating position, the width (W) of the product is the total distance from left to right; the depth (D) of the product is the total distance from the front to the back; the height (H) is the total distance from the bottom of the product to the top. Pictures may vary from actual product.

DISCLAIMER
- Output performance figures quoted are dependant on various factors. Scientific reserves the right, without notice, to make changes and revisions to product specifications, materials and design, when we believe it will provide better performance, durability, and efficiency. Unit measurements may vary depending on point of reference.

WARRANTY
- All Scientific equipment comes with a standard ONE year warranty, The 1 Year Extended Warranty is only valid for customers who have registered their warranty online at www.scientific.co.za within the first year after making their purchase.