
Scientific | Laboratory Equipment

PROUDLY MANUFACTURED IN SOUTH AFRICA

40L WATER BATH
Model: 132A



FEATURES

- Ideally suited to all general laboratory applications
- All baths together with the optional accessories are designed to be used with the Model 103 circulators
- Body is constructed from corrosion protected mild steel with an epoxy powder coated finish
- Inner tank is manufactured from grade 304 stainless steel for exceptional corrosion resistance
- 40L Capacity

NOTE

- Optional gabled lids and height reducers manufactured from grade 304 stainless steel

QUALITY GUARANTEED
ALL SCIENTIFIC EQUIPMENT COMES WITH A TWO YEAR WARRANTY ON COMPONENTS AND DEFECTIVE WORKMANSHIP - Ts & Cs APPLY.

Scientific | Laboratory Equipment

PROUDLY MANUFACTURED IN SOUTH AFRICA

40L WATER BATH
Model: 132A


SPECIFICATIONS

PACKAGING (W x D x H)	930 x 415 x 315 (mm)
UNIT (W x D x H)	835 x 350 x 245 (mm)
NET WEIGHT	15kg
GROSS WEIGHT	16kg
SHIPPING VOLUME	0.12m^3
TANK DIMENSIONS	705 x 305 x 200 (mm)
TANK VOLUME	40ℓ

OPTIONAL EXTRAS

CODE	DESCRIPTION
136A	Height Reducer Adjustable 50 - 90mm
139A	Stainless steel Gabled lid





NOTE
- When viewing the appliance from the front in it's usual operating position, the width (W) of the product is the total distance from left to right; the depth (D) of the product is the total distance from the front to the back; the height (H) is the total distance from the bottom of the product to the top. Pictures may vary from actual product.

DISCLAIMER
- Output performance figures quoted are dependant on various factors. Scientific reserves the right, without notice, to make changes and revisions to product specifications, materials and design, when we believe it will provide better performance, durability, and efficiency. Unit measurements may vary depending on point of reference.

WARRANTY
- All Scientific equipment comes with a standard ONE year warranty, The 1 Year Extended Warranty is only valid for customers who have registered their warranty online at www.scientific.co.za within the first year after making their purchase.

QUALITY GUARANTEED

ALL SCIENTIFIC EQUIPMENT COMES WITH A TWO YEAR WARRANTY ON COMPONENTS AND DEFECTIVE WORKMANSHIP - Ts & Cs APPLY.

www.scientific.co.za

Published Date: 2026.09.16