Scientific | Laboratory Equipment

MANUFACTURED IN SOUTH AFRICA

12.5mm RIFFLE SPLITTER

Model: 600



PROCESSING

FEATURES

- A sampling device used for sample splitting and acquiring a representative sample
- Battery of twenty open chambers of 12.5 mm
- Powder coated steel construction

NOTE

- Supplied with three bins
- Sample is equally divided into two subsamples

QUALITY GUARANTEED

ALL SCIENTIFIC EQUIPMENT COMES WITH A TWO YEAR WARRANTY ON COMPONENTS AND DEFECTIVE WORKMANSHIP - Ts & Cs APPLY.


Scientific | Laboratory Equipment

MANUFACTURED IN SOUTH AFRICA

12.5mm RIFFLE SPLITTER
Model: 600

SPECIFICATIONS

PACKAGING (W x D x H)	425 x 450 x 500 (mm)
UNIT (W x D x H)	425 x 450 x 450 (mm)
NET WEIGHT	8kg
GROSS WEIGHT	9kg
SHIPPING VOLUME	0.10m^3
CHAMBER SIZE	12.5 (mm)
NUMBER OF SLOTS	20



NOTE

- When viewing the appliance from the front in it's usual operating position, the width (W) of the product is the total distance from left to right; the depth (D) of the product is the total distance from the front to the back; the height (H) is the total distance from the bottom of the product to the top. Pictures may vary from actual product.

DISCLAIMER

- Output performance figures quoted are dependent on various factors. Scientific reserves the right, without notice, to make changes and revisions to product specifications, materials and design, when we believe it will provide better performance, durability, and efficiency. Unit measurements may vary depending on point of reference.

WARRANTY

- All Scientific equipment comes with a standard ONE year warranty, The 1 Year Extended Warranty is only valid for customers who have registered their warranty online at www.scientific.co.za within the first year after making their purchase.